



Promoting Excellence In Ultrasound

Policies and Statements

A3

Safety Statement On Continuous Wave Doppler Fetal Monitoring

Policies and Statements

A3

Safety Statement On Continuous Wave Doppler Fetal Monitoring

*September 1998, **Reaffirmed** February 2000, **Reaffirmed** July 2008*

Diagnostic devices emitting continuous wave Doppler ultrasound are used to monitor fetal and neonatal heartbeat.

This equipment emits low power levels and its use is not contra-indicated on safety grounds, even when applied for extended periods.



Australasian Society for Ultrasound in Medicine
PO Box 943, Crows Nest NSW 1585, Sydney AUSTRALIA
Phone: +61 2 9438 2078 • Fax: +61 2 9438 3686
Email: asum@asum.com.au • Website: <http://www.asum.com.au>
ABN: 64 001 679 161